DTB
BIOTECHNOLOGY

DTB BIOTECHNOLOGY LTD

Provides life time quarantee for implants



ABOUT DTB



Follow the light of technology

Founded by experienced dentists and expert researchers, DTB is conducting an intensive scientific R&D and design process in the fields of maxillofacial surgery, prosthetics, and periodontics to bring the UK's first 100% domestic implant brand to life.

Our development program is currently in the R&D, design, and prototyping phase. Our goal is to gain the trust of clinicians worldwide by bringing our scientifically validated, precision-engineered, and high-quality implant systems to life with 100% UK manufacturing.

DTB BIOTECHNOLOGY

Our Mission


Located at the heart of high technology in the United Kingdom, DTB's mission is to establish the country's independent and 100% domestic implant manufacturing ecosystem. With our high-precision product series currently in the design and R&D phase, we aim to attract patients traveling abroad for treatment back to local clinics and empower our clinics with top-quality domestic products.


What makes us different ?

Founded with the vision of being a premier 100% UK-based manufacturer, DTB BIOTECHNOLOGY LTD combines advanced R&D with scientific expertise to support modern medicine's philosophy of longevity.

While our core identity is a manufacturing enterprise, our projected product portfolio–dental implants, surgical kits, and biomaterials–is currently undergoing its final rigorous design, prototyping, and external biomechanical testing. Alongside our expert partners, we are meticulously preparing these pre-production systems for the manufacturing line. Driven by our goal of exclusive domestic production, we aim to eliminate foreign dependency and export high-quality UK biotechnology globally once mass production commences.

Expected PRODUCTS
SURFACE DESGIN





DTB Implants with an SLA Surface feature a combination of sand-blasting and acid-etching. This process increases the implant surface area by creating both micro and nano structures, allowing the surface to bind with more osteoblast cells. This surface also eliminates various surface contaminants through a thorough cleaning process, while acid-etching ensures a high-quality and pure implant surface.


The combination of sand-blasting and acid-etching has been proven to be successful in numerous clinical and scientific studies.

To ensure the highest standards of quality and objectivity, all advanced laboratory evaluations–including SEM, in-vitro, and sterilization studies–are rigorously conducted in collaboration with accredited external laboratories and expert business partners.

This process creates an osteoconductive morphology, achieving superb osseointegration results compared to other surface treatment methods. In addition to all these benefits, the DTB Active Surface maintains surface activation in its dry form as well, providing a fast osseointegration process and the shortest healing period.

Upcoming DTB Series

Pre-mass production prototype showcasing our precision engineering.



Micro Threads

Reduces bone stresses and minimize bone resorption

Taper Form

Cylindrical shape promotes long-term osseointegration by enlarging surface area and bone to implant contact

Self-Taping Threads

Ensure better primary stability and avoid excessive force causing cortical bone absorption and promotes bone tapering during insertion.

Internal Connection

Self-locking to enhance the connection between fixture and abutment with 11° Morse Taper Improves the stability and to prevent micro-movement. Prevents of micro leakage and bacteria intrusion. Perfect hex-connection for Ultra Mini, Mini.

Surface Treatment

SLA (Sand-Blasted and Acid-Etching) prosess for surface

Apical Part

Easy insertion
Helps prevent damage to anatomical structures

Upcoming DTB Series

Pre-mass production prototype showcasing our precision engineering.

Narrow and Bone **Level**

The DTB Biotechnology is incorporated in the Bone Level implant lines and ensures optimal esthetic results both in the anterior and posterior area.





Surgical **Kit**

The DTB Biotechnology is characterized by one unique surgical kit for bone level line. Its intuitive and elegant design prioritizes efficiency and ease of use.

Expected PRODUCTS
DTB IMPLANT SLA & ACT

AGRESSIVE SLA ULTRA MINI IMPLANT Ø3.0

Ref No	Size (L)	Ø	Material
DTBA3010UM	10,0 mm	Ø3.0	Titanium Grade 4
DTBA3011UM	11,5 mm	Ø3.0	Titanium Grade 4
DTBA3013UM	13,0 mm	Ø3.0	Titanium Grade 4
DTBA3015UM	15,0 mm	Ø3.0	Titanium Grade 4



AGRESSIVE SLA ACT ULTRA MINI IMPLANT Ø3.0

Ref No	Size (L)	Ø	Material
DTBA3010UMAct	10,0 mm	Ø3.0	Titanium Grade 4
DTBA3011UMAct	11,5 mm	Ø3.0	Titanium Grade 4
DTBA3013UMAct	13,0 mm	Ø3.0	Titanium Grade 4
DTBA3015UMAct	15,0 mm	Ø3.0	Titanium Grade 4



Expected PRODUCTS

DTB IMPLANT SLA & ACT

AGRESSIVE SLA MINI IMPLANT Ø3.5

Ref No	Size (L)	Ø	Material
DTBA3585M	8,5 mm	Ø3.5	Titanium Grade 4
DTBA3510M	10,0 mm	Ø3.5	Titanium Grade 4
DTBA3511M	11,5 mm	Ø3.5	Titanium Grade 4
DTBA3513M	13,0 mm	Ø3.5	Titanium Grade 4
DTBA3515M	15,0 mm	Ø3.5	Titanium Grade 4



AGRESSIVE SLA ACT MINI IMPLANT Ø3.5

Ref No	Size (L)	Ø	Material
DTBA3585MAct	8,5 mm	Ø3.5	Titanium Grade 4
DTBA3510MAct	10,0 mm	Ø3.5	Titanium Grade 4
DTBA3511MAct	11,5 mm	Ø3.5	Titanium Grade 4
DTBA3513MAct	13,0 mm	Ø3.5	Titanium Grade 4
DTBA3515MAct	15,0 mm	Ø3.5	Titanium Grade 4



Expected PRODUCTS
DTB IMPLANT SLA & ACT

AGRESSIVE SLA MINI IMPLANT Ø4.0

Ref No	Size (L)	Ø	Material
DTBA4070M	7,0 mm	Ø4.0	Titanium Grade 4
DTBA4085M	8,5 mm	Ø4.0	Titanium Grade 4
DTBA4010M	10,0 mm	Ø4.0	Titanium Grade 4
DTBA4011M	11,5 mm	Ø4.0	Titanium Grade 4
DTBA4013M	13,0 mm	Ø4.0	Titanium Grade 4
DTBA4015M	15,0 mm	Ø4.0	Titanium Grade 4



AGRESSIVE SLA ACT MINI IMPLANT Ø4.0

Ref No	Size (L)	Ø	Material
DTBA4070MAct	7,0 mm	Ø4.0	Titanium Grade 4
DTBA4085MAct	8,5 mm	Ø4.0	Titanium Grade 4
DTBA4010MAct	10,0 mm	Ø4.0	Titanium Grade 4
DTBA4011MAct	11,5 mm	Ø4.0	Titanium Grade 4
DTBA4013MAct	13,0 mm	Ø4.0	Titanium Grade 4
DTBA4015MAct	15,0 mm	Ø4.0	Titanium Grade 4



AGRESSIVE SLA MINI IMPLANT Ø4.5

Ref No	Size (L)	Ø	Material
DTBA4570M	7,0 mm	Ø4.5	Titanium Grade 4
DTBA4585M	8,5 mm	Ø4.5	Titanium Grade 4
DTBA4510M	10,0 mm	Ø4.5	Titanium Grade 4
DTBA4511M	11,5 mm	Ø4.5	Titanium Grade 4
DTBA4513M	13,0 mm	Ø4.5	Titanium Grade 4
DTBA4515M	15,0 mm	Ø4.5	Titanium Grade 4



AGRESSIVE SLA ACT MINI IMPLANT Ø4.5

Ref No	Size (L)	Ø	Material
DTBA4570MAct	7,0 mm	Ø4.5	Titanium Grade 4
DTBA4585MAct	8,5 mm	Ø4.5	Titanium Grade 4
DTBA4510MAct	10,0 mm	Ø4.5	Titanium Grade 4
DTBA4511MAct	11,5 mm	Ø4.5	Titanium Grade 4
DTBA4513MAct	13,0 mm	Ø4.5	Titanium Grade 4
DTBA4515MAct	15,0 mm	Ø4.5	Titanium Grade 4



Expected PRODUCTS
DTB IMPLANT SLA & ACT

AGRESSIVE SLA MINI IMPLANT Ø5.0

Ref No	Size (L)	Ø	Material
DTBA5070M	7,0 mm	Ø5.0	Titanium Grade 4
DTBA5085M	8,5 mm	Ø5.0	Titanium Grade 4
DTBA5010M	10,0 mm	Ø5.0	Titanium Grade 4
DTBA5011M	11,5 mm	Ø5.0	Titanium Grade 4
DTBA5013M	13,0 mm	Ø5.0	Titanium Grade 4
DTBA5015M	15,0 mm	Ø5.0	Titanium Grade 4



AGRESSIVE SLA ACT MINI IMPLANT Ø5.0

Ref No	Size (L)	Ø	Material
DTBA5070MAct	7,0 mm	Ø5.0	Titanium Grade 4
DTBA5085MAct	8,5 mm	Ø5.0	Titanium Grade 4
DTBA5010MAct	10,0 mm	Ø5.0	Titanium Grade 4
DTBA5011MAct	11,5 mm	Ø5.0	Titanium Grade 4
DTBA5013MAct	13,0 mm	Ø5.0	Titanium Grade 4
DTBA5015MAct	15,0 mm	Ø5.0	Titanium Grade 4



HEALING ABUTMENT ULTRA MINI

Ref No	Size (L)	Ø	Material
DTB2HLA3013UM	GH1 H3	Ø3.0	Titanium Grade 5
DTB2HLA3025UM	GH2 H5	Ø3.0	Titanium Grade 5
DTB2HLA3037UM	GH3 H7	Ø3.0	Titanium Grade 5



HEALING ABUTMENT MINI

Ref No	Size (L)	Ø	Material
DTB1HLA4513M	GH 1 H 3	Ø4.5	Titanium Grade 5
DTB1HLA4525M	GH 2 H 5	Ø4.5	Titanium Grade 5
DTB1HLA4537M	GH 3 H 7	Ø4.5	Titanium Grade 5
DTB1HLA5013M	GH 1 H 3	Ø5.0	Titanium Grade 5
DTB1HLA5025M	GH 2 H 5	Ø5.0	Titanium Grade 5
DTB1HLA5037M	GH 3 H 7	Ø5.0	Titanium Grade 5
DTB1HLA5513M	GH 1 H 3	Ø5.5	Titanium Grade 5
DTB1HLA5536M	GH 3 H 6	Ø5.5	Titanium Grade 5
DTB1HLA6013M	GH 1 H 3	Ø6.0	Titanium Grade 5
DTB1HLA6025M	GH 2 H 5	Ø6.0	Titanium Grade 5
DTB1HLA6037M	GH 3 H 7	Ø6.0	Titanium Grade 5









Expected PRODUCTS
DTB PROSTHETICS

SINUS HEALING ABUTMENT ULTRA MINI

Ref No	Size (L)	Ø	Material
DTBHLASWUM	–	Ø4.0	Titanium Grade 5



SINUS HEALING ABUTMENT MINI

Ref No	Size (L)	Ø	Material
DTBHLASWM	–	Ø4.5	Titanium Grade 5



ANATOMICAL STRAIGHT ABUTMENT ULTRA MINI

Ref No	Size (L)	Ø	Material
DTB2ESA301UM	GH:1,0 H:5,5	Ø3.0	Titanium Grade 5
DTB2ESA302UM	GH:2,0 H:5,5	Ø3.0	Titanium Grade 5
DTB2ESA303UM	GH:3,0 H:5,5	Ø3.0	Titanium Grade 5
DTB2ESA304UM	GH:4,0 H:5,5	Ø3.0	Titanium Grade 5
DTB2ESA305UM	GH:5,0 H:5,5	Ø3.0	Titanium Grade 5



STRAIGHT ABUTMENT MINI

Ref No	Size (L)	Ø	Material
DTB1SA45155M	GH:1,0 H:5,5	Ø4.5	Titanium Grade 5
DTB1SA45255M	GH:2,0 H:5,5	Ø4.5	Titanium Grade 5
DTB1SA45355M	GH:3,0 H:5,5	Ø4.5	Titanium Grade 5
DTB1SA45455M	GH:4,0 H:5,5	Ø4.5	Titanium Grade 5
DTB1SA45555M	GH:5,0 H:5,5	Ø4.5	Titanium Grade 5



STRAIGHT ABUTMENT MINI

Ref No	Size (L)	Ø	Material
DTB1SA50155M	GH 1 H 5.5	Ø5.0	Titanium Grade 5
DTB1SA50255M	GH 2 H 5.5	Ø5.0	Titanium Grade 5
DTB1SA50355M	GH 3 H 5.5	Ø5.0	Titanium Grade 5
DTB1SA50455M	GH 4 H 5.5	Ø5.0	Titanium Grade 5
DTB1SA50555M	GH 5 H 5.5	Ø5.0	Titanium Grade 5
DTB1SA55155M	GH 1 H 5.5	Ø5.5	Titanium Grade 5
DTB1SA55255M	GH 2 H 5.5	Ø5.5	Titanium Grade 5
DTB1SA55355M	GH 3 H 5.5	Ø5.5	Titanium Grade 5
DTB1SA55455M	GH 4 H 5.5	Ø5.5	Titanium Grade 5
DTB1SA55555M	GH 5 H 5.5	Ø5.5	Titanium Grade 5
DTB1SA55170M	GH 1 H 7.0	Ø5.5	Titanium Grade 5
DTB1SA55270M	GH 2 H 7.0	Ø5.5	Titanium Grade 5
DTB1SA55370M	GH 3 H 7.0	Ø5.5	Titanium Grade 5
DTB1SA55470M	GH 4 H 7.0	Ø5.5	Titanium Grade 5
DTB1SA55570M	GH 5 H 7.0	Ø5.5	Titanium Grade 5
DTB1SA60155M	GH 1 H 5.5	Ø6.0	Titanium Grade 5
DTB1SA60255M	GH 2 H 5.5	Ø6.0	Titanium Grade 5
DTB1SA60355M	GH 3 H 5.5	Ø6.0	Titanium Grade 5
DTB1SA60455M	GH 4 H 5.5	Ø6.0	Titanium Grade 5
DTB1SA60555M	GH 5 H 5.5	Ø6.0	Titanium Grade 5









Expected PRODUCTS
DTB PROSTHETICS

ANATOMICAL ANGLED ABUTMENT ULTRA MINI

Ref No	Size (L)	Ø	Material
DTB2ANA30215XUM	GH2 X-TYPE 15°	Ø3.0	Titanium Grade 5



ANATOMICAL ANGLED ABUTMENT MINI

Ref No	Size (L)	Ø	Material
DTB1EANA43215XM	GH2 X-TYPE 15°	Ø4.3	Titanium Grade 5
DTB1EANA43215YM	GH2 Y-TYPE 15°	Ø4.3	Titanium Grade 5
DTB1EANA43225XM	GH2 X-TYPE 25°	Ø4.3	Titanium Grade 5
DTB1EANA43225YM	GH2 Y-TYPE 25°	Ø4.3	Titanium Grade 5
DTB1EANA43415XM	GH4 X-TYPE 15°	Ø4.3	Titanium Grade 5
DTB1EANA43415YM	GH4 Y-TYPE 15°	Ø4.3	Titanium Grade 5
DTB1EANA43425XM	GH4 X-TYPE 25°	Ø4.3	Titanium Grade 5
DTB1EANA43425YM	GH4 Y-TYPE 25°	Ø4.3	Titanium Grade 5
DTB1EANA45215XM	GH2 X-TYPE 15°	Ø4.5	Titanium Grade 5
DTB1EANA45215YM	GH2 Y-TYPE 15°	Ø4.5	Titanium Grade 5
DTB1EANA45225XM	GH2 X-TYPE 25°	Ø4.5	Titanium Grade 5
DTB1EANA45225YM	GH2 Y-TYPE 25°	Ø4.5	Titanium Grade 5
DTB1EANA45415XM	GH4 X-TYPE 15°	Ø4.5	Titanium Grade 5
DTB1EANA45415YM	GH4 Y-TYPE 15°	Ø4.5	Titanium Grade 5
DTB1EANA45425XM	GH4 X-TYPE 25°	Ø4.5	Titanium Grade 5





ANATOMICAL ANGLED ABUTMENT MINI

Ref No	Size (L)	Ø	Material
DTB1EANA45425YM	GH4 Y-TYPE 25°	Ø4.5	Titanium Grade 5
DTB1EANA55215XM	GH2 X-TYPE 15°	Ø5.5	Titanium Grade 5
DTB1EANA55215YM	GH2 Y-TYPE 15°	Ø5.5	Titanium Grade 5
DTB1EANA55225XM	GH2 X-TYPE 25°	Ø5.5	Titanium Grade 5
DTB1EANA55225YM	GH2 Y-TYPE 25°	Ø5.5	Titanium Grade 5
DTB1EANA55415XM	GH4 X-TYPE 15°	Ø5.5	Titanium Grade 5
DTB1EANA55415YM	GH4 Y-TYPE 15°	Ø5.5	Titanium Grade 5
DTB1EANA55425XM	GH4 X-TYPE 25°	Ø5.5	Titanium Grade 5
DTB1EANA55425YM	GH4 Y-TYPE 25°	Ø5.5	Titanium Grade 5



BALL ATTACHMENT ULTRA MINI

Ref No	Size (L)	Ø	Material
DTB2BAS301UM	GH 1	Ø3.0	Titanium Grade 5
DTB2BAS302UM	GH 2	Ø3.0	Titanium Grade 5
DTB2BAS303UM	GH 3	Ø3.0	Titanium Grade 5
DTB2BAS304UM	GH 4	Ø3.0	Titanium Grade 5
DTB2BAS305UM	GH 5	Ø3.0	Titanium Grade 5



Expected PRODUCTS
DTB PROSTHETICS

BALL ATTACHMENT MINI

Ref No	Size (L)	Ø	Material
DTB1BAS351M	GH 1	Ø3.5	Titanium Grade 5
DTB1BAS352M	GH 2	Ø3.5	Titanium Grade 5
DTB1BAS353M	GH 3	Ø3.5	Titanium Grade 5
DTB1BAS354M	GH 4	Ø3.5	Titanium Grade 5
DTB1BAS355M	GH 5	Ø3.5	Titanium Grade 5



BALL ATTACHMENT ANALOG ULTRA MINI

Ref No	Size (L)	Ø	Material
DTB2BLAUM	–	–	Titanium Grade 5



BALL ATTACHMENT ANALOG MINI

Ref No	Size (L)	Ø	Material
DTB1BLAM	–	–	Titanium Grade 5



BALL ATTACHMENT CAP (NORMAL-HARD)

Ref No	Size (L)	Ø	Material
DTB1BACAPN	–	–	–
DTB1BACAPH	–	–	–



LOCATOR ABUTMENT ULTRA MINI

Ref No	Size (L)	Ø	Material
DTBLCT1UM	GH 1	–	Titanium Grade 5
DTBLCT2UM	GH 2	–	Titanium Grade 5
DTBLCT3UM	GH 3	–	Titanium Grade 5
DTBLCT4UM	GH 4	–	Titanium Grade 5
DTBLCT5UM	GH 5	–	Titanium Grade 5



LOCATOR ABUTMENT MINI

Ref No	Size (L)	Ø	Material
DTBLCT1M	GH1	–	Titanium Grade 5
DTBLCT2M	GH2	–	Titanium Grade 5
DTBLCT3M	GH3	–	Titanium Grade 5
DTBLCT4M	GH4	–	Titanium Grade 5
DTBLCT5M	GH5	–	Titanium Grade 5



LOCATOR ABUTMENT ANALOG

Ref No	Size (L)	Ø	Material
DTBLILA	–	–	Titanium Grade 5



LOCATOR ABUTMENT TRANSFER COPING

Ref No	Size (L)	Ø	Material
DTB1LOCITI	–	–	Titanium Grade 5



Expected PRODUCTS
DTB PROSTHETICS

LOCATOR ABUTMENT CAP (NORMAL -HARD)

Ref No	Size (L)	Ø	Material
DTB1LCAPN	–	–	–
DTB1LCAPH	–	–	–



IMPLANT LAB. ANALOG ULTRA MINI

Ref No	Size (L)	Ø	Material
DTB2ILAUM	–	–	Titanium Grade 5



IMPLANT LAB. ANALOG MINI

Ref No	Size (L)	Ø	Material
DTB1ILAM	–	–	Titanium Grade 5



DIGITAL ANALOG ULTRA MINI

Ref No	Size (L)	Ø	Material
DTBDGAUM	–	–	Titanium Grade 5




DIGITAL ANALOG MINI

Ref No	Size (L)	Ø	Material
DTBDGAM	–	–	Titanium Grade 5



OPEN TRAY TRANSFER COPING ULTRA MINI

Ref No	Size (L)	Ø	Material
DTB1AKITIUMH (Hex)	–	–	Titanium Grade 5
DTB1AKITIUM (Non-Hex)	–	–	Titanium Grade 5



OPEN TRAY TRANSFER COPING MINI

Ref No	Size (L)	Ø	Material
DTB1AKITIMH (Hex)	–	–	Titanium Grade 5
DTB1AKITIM (Non-Hex)	–	–	Titanium Grade 5



CLOSED TRAY TRANSFER COPING ULTRA MINI

Ref No	Size (L)	Ø	Material
DTB2ITI3011UM	11	Ø3.0	Titanium Grade 5
DTB2ITI3014UM	14	Ø3.0	Titanium Grade 5



CLOSED TRAY TRANSFER COPING MINI

Ref No	Size (L)	Ø	Material
DTB1ITI4011M	11	Ø4.0	Titanium Grade 5
DTB1ITI4014M	14	Ø4.0	Titanium Grade 5
DTB1ITI4511M	11	Ø4.5	Titanium Grade 5
DTB1ITI4514M	14	Ø4.5	Titanium Grade 5
DTB1ITI5011M	11	Ø5.0	Titanium Grade 5
DTB1ITI5014M	14	Ø5.0	Titanium Grade 5



Expected PRODUCTS
DTB PROSTHETICS

PEEK ABUTMENT ULTRA MINI

Ref No	Size (L)	Ø	Material
DTB1PAPUMH (Hex)	–	–	PEEK
DTB1PAUM (Non-Hex)	–	–	PEEK



PEEK ABUTMENT MINI

Ref No	Size (L)	Ø	Material
DTB1PAMH (Hex)	–	–	PEEK
DTB1PAM (Non-Hex)	–	–	PEEK



UNIVERSAL ABUTMENT ULTRA MINI

Ref No	Size (L)	Ø	Material
DTB1UAUMH (Hex)	–	–	Titanium Grade 5
DTB1UAUM (Non-Hex)	–	–	Titanium Grade 5



UNIVERSAL ABUTMENT MINI

Ref No	Size (L)	Ø	Material
DTB1UAMH (Hex)	–	–	Titanium Grade 5
DTB1UAM (Non-Hex)	–	–	Titanium Grade 5



Ti-BASE ABUTMENT ULTRA MINI

	Ref No	Size (L)	Ø	Material
(Hex)	DTBTBASE1UMHS	GH 0.6 S	–	Titanium Grade 5
	DTBTBASE2UMHS	GH 1.5 S	–	Titanium Grade 5
	DTBTBASE3UMHS	GH 2.5 S	–	Titanium Grade 5
(Non-Hex)	DTBTBASE1UMS	GH 0.6 S	–	Titanium Grade 5
	DTBTBASE2UMS	GH 1.5 S	–	Titanium Grade 5
	DTBTBASE3UMS	GH 2.5 S	–	Titanium Grade 5





Ti-BASE ABUTMENT MINI

	Ref No	Size (L)	Ø	Material
(Hex)	DTBTBASE1MHS	GH 0.6 S	–	Titanium Grade 5
	DTBTBASE2MHS	GH 1.5 S	–	Titanium Grade 5
	DTBTBASE3MHS	GH 2.5 S	–	Titanium Grade 5
(Non-Hex)	DTBTBASE1MS	GH 0.6 S	–	Titanium Grade 5
	DTBTBASE2MS	GH 1.5 S	–	Titanium Grade 5
	DTBTBASE3MS	GH 2.5 S	–	Titanium Grade 5





Ti-BASE ABUTMENT CASTABLE MANUAL

Ref No	Size (L)	Ø	Material
DTBTCASM	–	–	Titanium Grade 5



PREMILL ABUTMENT ULTRA MINI

Ref No	Size (L)	Ø	Material
DTBPREUMA	–	–	Titanium Grade 5



Expected PRODUCTS
DTB PROSTHETICS

PREMILL ABUTMENT MINI

Ref No	Size (L)	Ø	Material
DTBPREMA	–	–	Titanium Grade 5



TELESCOPIC ABUTMENT ULTRA MINI

Ref No	Size (L)	Ø	Material
DTB1TELAUM	–	–	Titanium Grade 5



TELESCOPIC ABUTMENT MINI

Ref No	Size (L)	Ø	Material
DTB1TELAM	–	–	Titanium Grade 5



UCLA ABUTMENT MINI

Ref No	Size (L)	Ø	Material
DTB1UCLAMH (Hex)	–	–	Pom C-Delrin
DTB1UCLAM (Non-Hex)	–	–	Pom C-Delrin



MULTI UNIT STRAIGHT ABUTMENT MINI

Ref No	Size (L)	Ø	Material
DTB1MUSA1M	GH1	–	Titanium Grade 5
DTB1MUSA3M	GH3	–	Titanium Grade 5
DTB1MUSA5M	GH5	–	Titanium Grade 5



MULTI UNIT ANGLED ABUTMENT MINI

Ref No	Size (L)	Ø	Material
DTB1MUANA201M	GH1 20°	–	Titanium Grade 5
DTB1MUANA203M	GH3 20°	–	Titanium Grade 5
DTB1MUANA205M	GH5 20°	–	Titanium Grade 5
DTB1MUANA351M	GH1 35°	–	Titanium Grade 5
DTB1MUANA353M	GH3 35°	–	Titanium Grade 5
DTB1MUANA355M	GH5 35°	–	Titanium Grade 5





MULTI UNIT HEALING ABUTMENT MINI

Ref No	Size (L)	Ø	Material
DTB1MUHLA	–	–	Titanium Grade 5



MULTI UNIT OPEN TRAY TRANSFER COPING MINI

Ref No	Size (L)	Ø	Material
DTB1MUITIH (Hex)	–	–	Titanium Grade 5
DTB1MUITI (Non-Hex)	–	–	Titanium Grade 5



Expected PRODUCTS
DTB PROSTHETICS

MULTI UNIT CLOSED TRAY TRANSFER COPING MINI

Ref No	Size (L)	Ø	Material
DTB1MUITIC	–	–	Titanium Grade 5



TEMPORARY ABUTMENT MINI

Ref No	Size (L)	Ø	Material
DTB1MUTAH (Hex)	–	–	Titanium Grade 5
DTB1MUTA (Non-Hex)	–	–	Titanium Grade 5



MULTI UNIT CASTABLE ABUTMENT MINI

Ref No	Size (L)	Ø	Material
DTB1MUDAH (Hex)	–	–	Pom C-Delrin
DTB1MUDA (Non-Hex)	–	–	Pom C-Delrin



MULTI UNIT Ti-BASE ABUTMENT MINI

Ref No	Size (L)	Ø	Material
DTBMTBASE	–	–	Titanium Grade 5



MULTI UNIT STRAIGHT ANALOG MINI

Ref No	Size (L)	Ø	Material
DTB1MUILAH (Hex)	–	–	Titanium Grade 5
DTB1MUILA (Non-Hex)	–	–	Titanium Grade 5



MULTI UNIT ANGLED 20° - 35° ANALOG MINI

Ref No	Size (L)	Ø	Material
DTB1MUILAA20	–	–	Titanium Grade 5
DTB1MUILAA35	–	–	Titanium Grade 5



MULTI DIGITAL ANALOG MINI

Ref No	Size (L)	Ø	Material
DTBDGAMULTI	–	–	Titanium Grade 5



MULTI UNIT CARRIER PART MINI

Ref No	Size (L)	Ø	Material
DTB1MUYP	–	–	Titanium Grade 5



Expected PRODUCTS
DTB PROSTHETICS

MULTI UNIT STRAIGHT ANALOG MINI

Ref No	Size (L)	Ø	Material
DTBMSCP	–	–	Titanium Grade 5



EXOCAD SCAN BODY

Ref No	Size (L)	Ø	Material
DTBESCANL (Large)	–	–	Pom C-Delrin
DTBESCANS (Small)	–	–	Pom C-Delrin



IMPLANT SCAN BODY

Ref No	Size (L)	Ø	Material
DTBSCANPUM	–	–	Titanium Grade 5
DTBSCANPM	–	–	Titanium Grade 5



CEREC SCAN BODY

Ref No	Size (L)	Ø	Material
DTBTSCANL (Large)	–	–	ABS
DTBTSCANS (Small)	–	–	ABS




Expected PRODUCTS
SURGICAL KIT



DENTAL DRILL


Ref No	Ø	Material
DTBDLADR (Lance Drill)	–	Stainless Steel
DTBFAG23	Ø2.3 AG	Stainless Steel
DTBFAG30	Ø3.0 AG	Stainless Steel
DTBFAG35	Ø3.5 AG	Stainless Steel
DTBFAG40	Ø4.0 AG	Stainless Steel
DTBFAG45	Ø4.5 AG	Stainless Steel
DTBFAG50	Ø5.0 AG	Stainless Steel



Expected PRODUCTS
DTB SURGICAL KIT

CORTICAL DENTAL DRILL

Ref No	Ø	Material
DTBFAG30C	Ø3.0 AG	Stainless Steel
DTBFAG35C	Ø3.5 AG	Stainless Steel
DTBFAG40C	Ø4.0 AG	Stainless Steel
DTBFAG45C	Ø4.5 AG	Stainless Steel
DTBFAG50C	Ø5.0 AG	Stainless Steel



DRIVERS

Ref No	Ø	Material
Driver 1.2 S	DTBDR12S	Stainless Steel
Driver 1.2 L	DTBDR12L	Stainless Steel
Driver Torque 1.2 S	DTBDT12S	Stainless Steel
Driver Torque 1.2 L	DTBDR12L	Stainless Steel



EXTENSIONS


Name	Ref No	Material
Mount Extension S	DTBDMEXS	Stainless Steel
Mount Extension L	DTBDMEXL	Stainless Steel
Drill Extension	DTBDDREX	Stainless Steel



Expected PRODUCTS
DTB SURGICAL KIT

MOUNT DRIVERS

Name	Ref No	Material
Mount Driver S	DTBDMODS	Stainless Steel
Mount Driver L	DTBDMODL	Stainless Steel



PARALEL PIN

Name	Ref No	Material
Parallel Pin	DTBDPLP4	Stainless Steel
Parallel Pin Angled 20°	DTBDPLP20	Stainless Steel
Parallel Pin Angled 35°	DTBDPLP35	Stainless Steel




REMOVAL TOOLS

Name	Ref No	Material
Removal Tool Ultra Mini	DTBDRTOUM	Stainless Steel
Removal Tool Mini	DTBDRTOM	Stainless Steel



FIXTURE DRIVERS

Name	Ref No	Material
Fixture Driver Ultra Mini	DTBDFDRUM	Stainless Steel
Fixture Driver Mini S	DTBDFDRMS	Stainless Steel
Fixture Driver Mini L	DTBDFDRML	Stainless Steel
Depht Gauge	DTBDEPTH	Stainless Steel
Ratchet	DTBDRACHT	Stainless Steel



Bovine Collagen Membrane

DTB Bovine Collagen Membrane is extracted from the bovine pericardium and is sterilized using gamma radiation. It is used for the strengthening, regeneration and replacement of connective tissues.
DTB Bovine Collagen Membrane initially covers the tissue defect safely and acts as a template for the cellular repair mechanism that replaces the structure with the body's own natural tissues.



Membrane		
15x20 mm	20x30 mm	30x40 mm

Bovine Bone Graft

DTB Bovine Bone Graft is obtained from the bovine cancellous bone and is subjected to gamma sterilisation.
DTB Bovine Bone Graft is used to fill bone defects, for establishing stabilisation and for plastic reconstruction of an extracted or damaged bone section.



Cancellous - Cortical			
0.5 cc	1 cc	2 cc	5 cc
0.5 cc	0.25 - 1 mm	0.25 - 1 mm	0.25 - 1 mm
1 - 2 mm	1 - 2 mm	1 - 2 mm	1 - 2 mm

*Halal certified.

DTB
BIOTECHNOLOGY

DTB BIOTECHNOLOGY LTD

Provides life time quarantee for implants



DTB
BIOTECHNOLOGY

DTB BIOTECHNOLOGY LTD

Provides life time guarantee for implants

DTB BIOTECHNOLOGY LTD
Company Number: 17062348
Unit A3, Gateway Tower, 32 Western Gateway, London, United Kingdom, E16 1YL

dtbbiotechnology.com